

ASPROMUS[®]

Innovative Medical Aesthetic Equipment

2404B

Medical Aesthetic Treatment Chair

Designed for injectables and advanced aesthetics



CE FDA RoHS

Aspromus, a division of Silverfox

www.aspromus.us

info@silverfoxamerica.com

The Aspromus Difference

Thoughtfully engineered for injectables, lasers, and precision-driven clinical workflows.

Upright Injection Position (Up to 82°)

Supports steady, comfortable patient alignment for precise injectable work. Ideal for neuromodulators, dermal fillers, and high-precision facial treatments.



Sleekly Slim base. Securely Grounded

Experience the perfect fusion of form and function. Our engineered ultra-slim base offers a sophisticated silhouette while allowing practitioners to work close without obstruction or toe-bumping. Beneath its elegant profile lies a robust structure ensuring absolute stability, even during full motorized adjustment.



Streamlined Control System

Side-mounted controls make it easy to adjust the backrest and height without visual clutter. A clean, simplified setup minimizes clutter, allowing the provider to stay focused on the procedure and the patient experience.



Medical-Grade, Detachable Upholstery

Silicone leather resists staining from CHG prep, numbing creams, and marking pens. Cushions detach for easy cleaning. Supports strict hygiene requirements and fast patient turnover.



Intelligent Memory System

Memory presets and reset function support repeatable injector and aesthetician ergonomics. Saves time between treatments and maintains a consistent posture.



Designed for the Demands of Aesthetic Medicine




Trendelenburg Position (When Needed)

Useful for managing vasovagal episodes and select intimate wellness treatments. Provides added clinical versatility and peace of mind, especially in practices that handle occasional fainting episodes or offer pelvic and rejuvenation services.

Specifications

Voltage :	AC 100-120V or AC 220-240V
Frequency :	50 / 60Hz
Max load bearing :	200 kg / 440 lbs
Net weight:	110 kg / 243 lbs
Gross weight:	123 kg / 273 lbs
Package size:	1400 x 810 x 610 (mm)≈0.69 m ³ 55 x 32 x 24 (in)≈24.5 cu ft
Product size(L*W*H) :	1810 x 859 x 570(min)/820(max) mm 71 x 34 x 22.4(min)/33(max) in

Standard Configuration





Options



Unit

mm / in



➤ Premium Silicone Upholstery

Pebble Grain



➤ Medical-grade, Stain-resistant Silicon Leather

Our silicone leather combines the premium feel of real leather with unmatched, easy-clean durability.

Effortless Cleanup: Simply wipe away everyday spills with mild soap and water. That's all it takes.

Bleach-Safe: Disinfect thoroughly with a diluted bleach solution (10%) to eliminate germs.

Alcohol-Resistant: Wipe down with 75% alcohol sanitizers for quick cleaning without damage.

➤ Sponge Dissection Diagram

Leather

Low Hardness Foam

High Firmness Foam

Laminate Support



*Actual colors may differ slightly from those displayed in images.

➤ Specifications

•Material Composition

Surface: 100% Silicone
Backing: Polyester

•Inherent Flame Resistance

FMVSS 302
BS5852 Crib 5

•Color Fastness

Light AATCC16.3 200h: Grade 3-4
Crocking AATCC8 Wet & Dry: Grade 5
Perspiration AATCC15: Grade 5

•Stain Resistance

CFFA-141,Excellent

•Physical Performance

Wyzenbeek Abrasion:100,000 double rubs
UV Stability ASTM D4329-05
Flex Resistance ASTM D2097-91
Hydrolysis Resistance ASTM D3690-02
Cold Crack Resistance CFFA-6
Salt Spray Resistance ASTM B117
Low-temp Flex Resistance ISO 17694
Odorless

•Sustainability

VOCs Free
Recyclable
RoHS 2.0 / REACH / Prop 65
Advanced Solvent Free Technology
No Toxic By-products
No Formaldehyde, Lead, PFC,Antimony Trioxide, Tri
No Plasticizers, Phthalates,Heavy Metals, BPA
No Harm to Air or Water

•Biocompatibility

Repeated Stimulation: ISO 10993-10:2010 Level 0
Sensitization: ISO 10993-10:2010 Level 0
Cytotoxicity: ISO 10993-5:2009 Level 1

